









Alois

Maximal comfort and
a minimal height



Technical information Alois

							
width	length	position min.	position max.	inclination head-support	inclination knee-buckle	inclination leg-support	users weight
107 cm	207 cm	20 cm	72 cm	65°	16°	30°	140 kg

Deviations in sizes and weight (1,5 cm, 1,5 kg, 1,5°) are possible



Comfortable and practical

The electronic controller easy allows you to adjust the:

- Height
- Head support
- Knee buckle (the foot part can be adjusted manually)

All 4 wheels can be blocked separately. The bed can easily be disassembled and transported. Side panels and self-raising device (4 possible positions) are part of the standard equipment.

Adaptable to your needs

For the past 50 years, Vermeiren has played an important role in home care products. We continually strive to meet your individual needs. Our products are the result of continuous research and many years of experience.

N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

Visit our website: www.vermeiren.com



VERMEIREN
Opening up new horizons in mobility

For more information, contact your local dealer

