


Detail A




Material			
No.	Plan	Article	Number
1	X004389	1903684 (Small)	1
	X003883	1901986 (Medium)	1
	X003884	1901987 (Large)	1
	P00079637	1991264 (XXL)	1
2	Y002759	1217101 (Large)	1
	Y003430	1217001 (Small)	1
3	Y011781	1678400	1
4	Y011782	1678501	1
5	X001101	1001605	2
6	---	1000973	4
7	---	1001564	4
8	---	1001589	2
9	---	1003021	4
10	---	1006999	4

Seat width	Combination
39	A
42	B
44	B
46	C
48	D
50	D
55	D
60*	E

- A: 1903684 + 1217101
- B: 1901986 + 1217001
- C: 1901986 + 1217101
- D: 1901987 + 1217101
- | E: 1991264 + 1217101

* - use with EZ1 size:57

B12 ALL MODELS				Plan No.
				E000538/D
Scale: 1:4	Drafted	Verified	Confirmed	
Date	18/12/2008	12/04/2024	12/04/2024	
Init.	SD	AZ	AS	