



Material			
No.	Plan	Article	Quantity
1	XP02831 (High 100)	2994640	1
	XP02832 (High 120)	2994641	1
2	XP02833 (High 100)	2994642	1
	XP02834 (High 120)	2994643	1
3	SANMED C-460/7000 Grey (Vinterno)	2993898	5,0m
	SANMED C-460/5202 Blue (Vermeiren)	2994651	5,0m
4*	Slide fastener (Vinterno)	2994224	5,0m
	Slide fastener (Vermeiren)	1995604	5,0m
5*	Slide fastener head (Vinterno)	2994225	2
	Slide fastener head (Vermeiren)	1995605	2
6*	Sewing thread (Vinterno)	W200319	20m
	Sewing thread (Vermeiren)	W200016	20m

Final dimensions of mattress: 1900x850x100mm or 1900x850x120mm
 Final dimensions of mattress extra: 235x850x100mm or 235x850x120mm /
 Wymiary końcowe materaca: 1900x850x100mm lub 1900x850x120mm
 Wymiary końcowe dostawki materaca: 235x850x100mm lub 235x850x120mm

EP00467 / A

* not depicted / nie ujęto na rysunku

MATTRESS FOR ILLICO BEDS				Plan No.
				EP00467 / A
Scale: 1:15	Drafted	Verified	Confirmed	
Date	28/10/2019	04/02/2021	04/02/2021	
Initials	AZ	EB	TK	